

Monnal T50

Breathing naturally



monnal T50

		Pressure Modes				Volume Modes			
		PSV		APCV		AVCV		SIMV	
		Adult	Child	Adult	Child	Adult	Child	Adult	Child
Settings		Invasive Ventilation / Non Invasive Ventilation Simple & Dual Limb							
Volume	VT (mL)					100 - 2000	50 - 500	100 - 2000	50 - 500
	Flow pattern	constant or decelerated							
Pressure	PI (cm H2O)	5 - 50		5 - 50				5 - 50	
	Slope	1 - 4		1 - 4				1 - 4	
	PEEP (cm H2O)	0 - 20							
Time	RRmin (bpm)	5 - 40	5 - 60	5 - 40	5 - 60	5 - 40	5 - 60		
	RRsimv (bpm)							1 - 40	1 - 60
	I:E					1 - 9		1 - 9	
	Ti (s)							0.3 - 3	
Triggers	Ti max (s)	0.3 - 3						0.3 - 3	
	Trig I (L/min)	0.5 - 10		(off) 0.5 - 10		(off) 0.5 - 10		0.5 - 10	
Sigh	Trig E (%)	10 - 90						10 - 90	
	Sigh period			9 - 200		9 - 200			
Apnea	Sigh coeff			1 - 2		1 - 2			
	T Apnea (s)	10 - 60	4 - 60					10 - 60	4 - 60
	RR Apnea (bpm)	7 - 40	9 - 60					7 - 40	9 - 60
	Vt Apnea (mL)	100 - 2000	50 - 500					100 - 2000	50 - 500

Specifications

Weight	5.3 kg
Size	33 x 25 x 18 mm
Screen	6"4
USB	2 slots

Specific modes

Programs	2
VT Target	✓
Boost Function	✓

Monitoring

Trends	48h
---------------	-----

Power supply

External AC	100 - 240 V, 50 to 60 Hz, 100 VA
External DC	10 - 26 V, 10 A Max
Internal battery	4h to 6h*
External battery	8h to 12h*

* Depending on settings



Monnal T50
KC027500



Carrying Bag
KF007700



Roll Stand
KC027700



External Battery
KC027400



Expiratory Valve SUx1
KY677400 + KY632500

For all references, please contact us



Air Liquide's Healthcare business supplies medical gases & equipment to hospitals, medical services to patients at home and is active in the fight against nosocomial infections. Its 8,000 employees serve 6,000 hospitals and 500,000 patients, on five continents. Air Liquide is the world leader in gases for industry, health and the environment.

Parc de Haute Technologie - 6, rue Georges Besse - 92182 Antony Cedex - Tél: +33 (0)1 40 96 66 00 - Fax : +33 (0)1 40 96 67 00